

PROJECT NAME : ProRenew Complex CLR



FORMULA NAME : Regenerating Serum

TRIAL CODE : 421.002.01

Lab Formula INCI Supplier

<i>Phase</i>	<i>Material Name</i>	<i>Mixed INCI</i>	<i>Supplier</i>	<i>% Material</i>
A	DEIONIZED WATER	WATER	N.A.	90.55
A	PEMULEN TR-1	ACRYLATES/C10-30 ALKYL ACRYLATE CROSSPOLYMER	LUBRIZOL	0.20
A	VERSTATIL TBG	TRIETHYLCITRATE GLYCERYL CAPRYLATE BENZOIC ACID	DR STRAETMANS	1.00
B	PHYTODEFENSE CLR	GLYCINE SOJA (SOYBEAN) OIL DICAPRYLYL ETHER MAGNOLIA GRANDIFLORA BARK EXTRACT LAURYL ALCOHOL	CLR	3.00
B	VITAMIN F FORTE	LINOLEIC ACID LINOLENIC ACID	CLR	1.00
B	LEXFEEL NATURAL	HEPTYL UNDECYLENATE	INOLEX	1.00
B	COPHEROL F 1300 C	TOCOPHEROL	BASF	0.05
C	PRORENEW COMPLEX CLR	LACTOCOCCUS FERMENT LYSATE	CLR	3.00
D	NAOH 32%	WATER SODIUM HYDRIXIDE	N.A	0.20
				100.00

Operating instructions

Operating Instructions

Mix A and stir slowly until completely soaked. Add B and homogenize with Ultra Turrax for 2 minutes. Add C while stirring. Adjust pH value to 5 - 5.5 with D.

The recommendations and formulations given are based on our knowledge and experience in the field of technical application.

They are, to the best of our belief, correct, but are offered without obligation.

Those who use our recommendations and formulations as well those as those who process CLR Active Agents are themselves responsible for the adherence to prevailing statutory regulations and the observance of patent rights as well as other protective rights for other companies.

This formula has been manufactured and stability tested using a special preservative, but has not been subjected to microbiological challenge tests.