

Moisturizing Hair Mask (O/W)



DayMoist CLR™

FORMULA No. 148.07.0877

Phase	Trade Name	INCI Name	w/w %	Supplier
A	Water	Water	ad 100	
B	Lanette O	Cetearyl Alcohol	4.00	BASF
	Kerabase	Brassica Alcohol, Brassicamidopropyl Dimethylamine, Polyester-11, Aspartic Acid, EDTA	15.00	Inolex
C	Microcare SB*	Water, Sodium Benzoate, Potassium Sorbate	1.00	Thor GmbH
	DayMoist CLR™	Water, Hydrolyzed Corn Starch, Beta Vulgaris (Beet) Root Extract	3.00	CLR

Manufacture:

Heat up A and B to 80°C-85°C separately. Add A to B and cool down while stirring. Add C in the given order at room temperature and mix until homogenous.

The recommendations and formulations given are based on our knowledge and experience in the field of technical application. The are, to the best of our belief, correct, but are offered without obligation. Those who use our recommendations and formulations as well as those who process CLR Active Agents are themselves responsible for the adherence to prevailing statutory regulations and the observance of patent rights as well as other protective rights for other companies.

*) This formula has been manufactured and stability tested using a special preservative, but has not been subjected to microbiological challenge tests.