



Tanning Drops
II-48.385.03

Phase	Material Name	US INCI	Supplier	% Material
A	Deionized Water	Water	N.A.	63.9
A	Glycerol 85%, vegetable Ph. Eur.	Glycerin	Gustav Heess	5.0
A	Keltrol T Plus	Xanthan Gum	CP Kelco	0.2
A	MultiMoist CLR™	Fructooligosaccharides Beta Vulgaris (Beet) Root Extract Water	CLR	3.0
A	Euxyl K 712	Sodium Benzoate Potassium Sorbate Water	Ashland	0.8
A	Citric Acid (50%)	Citric Acid Water	N.A.	0.6
B	Deionized Water	Water	N.A.	15.0
B	Dihydroxyacetone extra pure	Dihydroxyacetone	Merck	7.5
B	Erythrulose	Erythrulose Water	DSM	4.0
				100.0

Operating Instructions

Mix A and stir until uniform. Premix B and add it to A while stirring. Adjust pH value to 4, if necessary.

Directions for use:

Apply to clean, dry skin that has been exfoliated. Mix 2-3 tanning drops with your day or night moisturizer in the palm of your hand. Avoid eyebrows and hairline. Wash hands with soap and water immediately after use.

The recommendations and formulations given are based on our knowledge and experience in the field of technical application.

They are, to the best of our belief, correct, but are offered without obligation.

Those who use our recommendations and formulations as well as those who process CLR Active Agents are themselves responsible for the adherence to prevailing statutory regulations and the observance of patent rights as well as other protective rights for other companies.

This formula has been manufactured and stability-tested using a special preservative, but has not been subjected to microbiological challenge tests.

CLR - Chemisches Laboratorium Dr. Kurt Richter GmbH - www.clr-berlin.com